


SMMS - PE/PP NON-WOVEN UROLOGY GOWN



Recommended for all urological procedures and any other procedure with a large quantity of fluids.

- front and sleeves in fully waterproof PE/PP non-woven fabric
- back in breathable water-repellent SMMS non-woven fabric
- model identification collar: RED
- optimal coverage of operator in sitting position
- 4 internal/external ties
- ultrasonic sealed
- 2 hand towels
- **wrapped in sterile wrap. Dispenser box - Sterile**

CODE	SIZE size	LENGHT cm	STERILE	SECONDARY PACKAGING pieces per box
503008AHBSME	M	120	YES	30
503009AHBSME	L	130	YES	30
503010AHBSME	XL	140	YES	30
503011AHBSME	XXL	150	YES	30

NON-STERILE SURGICAL GOWNS

LIGHT SMMS NON-WOVEN "SP" SPECIALIST GOWN



- with 2 hand towels
- wrapped in sterile wrap
- **non sterile**

CODE	SIZE size	LENGHT cm	STERILE	SECONDARY PACKAGING pieces per box
500013CC3NMD	M	115	NO	50
500014CC3NMD	L	130	NO	50
500015CC3NMD	XL	135	NO	50
500016CC3NMD	XXL	145	NO	50

LIGHT SMMS NON-WOVEN "SP" REINFORCED GOWN



- with 2 hand towels
- wrapped in sterile wrap
- **non sterile**

CODE	SIZE size	LENGHT cm	STERILE	SECONDARY PACKAGING pieces per box
501013CC3NMD	M	115	YES	50
501014CC3NMD	L	130	YES	50
501015CC3NMD	XL	135	YES	50
501016CC3NMD	XXL	145	YES	50